



Test Certificate

Part Basket Holder with Stirring Shaft
Part Number 07558-06
Serial Number N18002992 - N18003628



| | | | USP / Ph.Eur. Tolerances |
|-------------------|---|---------------|-----------------------------|
| Dimensions | Shaft Ø | 9.4 – 10.1 mm | - - |
| | Cap Ø | 25.0 mm | ± 3.0 mm |
| | Cap size | 5.1 mm | ± 0.5 mm |
| | Clear opening | 20.2 mm | ± 1.0 mm |
| | Vent hole Ø | 2.0 mm | ± 0.5 mm |
| Material | Stainless Steel (electrolytic polished) | | |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2019-02-22

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland