




Part Vessel Mini Stirrer
 Part Number 10583
 Serial Number 78946 - 78966



		Tolerances		
Dimensions	Vessel	Outside Ø A	124.5 mm	+ 0.3 / - 0.2 mm
		Inside Ø B	100.0 mm	± 0.5 mm
		Outside Ø C	50.0 mm	± 1.0 mm
		Inside Ø D	40.0 mm	± 1.0 mm
		Height E	211.5 mm	± 0.5 mm
		With hemispherical bottom		
Material	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

Aesch, 2019-05-16

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate