




Part 1 Liter Vessel
 Part Number 14143-11
 Serial Number P35 001101 – P35 001062



**USP / Ph.Eur.
Tolerances**

Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	98 – 106 mm
		Height b With hemispherical bottom	168.2 mm	160 – 210 mm
Material	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2021-05-05

05-05-2021 14:51:16

SOTAX Quality Management
 EL 142/143 | Ground Floor | TTC Industrial area | MIDC Mahape | Navi Mumbai 400701 | India

Test Certificate