




Test Certificate

Part 1 Liter Vessel Inactinic Brown
Part Number 14143-12
Serial Number P35 000135 – P35 000190



		USP / Ph.Eur. Tolerances		
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	98 – 106 mm
		Height b	168.2 mm	160 – 210 mm
		With hemispherical bottom		
Material	Borosilicate glass inactinic brown			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2022-04-25

25-04-2022 15:02:34

Moreshwar Pulekar
SOTAX Quality Management
C-64 | First Floor | TTC Industrial area | Turbhe | Navi Mumbai 40070 | India