


Test Certificate

Part: **China Vessel 250 mL**
 Number: 14536-01
 Serial Number: 66720 - 66751

For SOTAX AT Dissolution Testing Systems



Dimensions:

| | | |
|----------|---|--------|
| Inside Ø | A | 62 mm |
| Height | B | 126 mm |

PPRC
Tolerances
± 3 mm
± 6 mm

With hemispherical bottom

Material:

Borosilicate glass

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with PPRC Pharmacopoeia requirements (at time of delivery).

SOTAX Quality Management
Bruno Balli
2016-05-12


Signature