

# AT7 *Smart*

## AT 7*smart* On-line Automated Dissolution System

SWISS  
QUALITY  
ISO 9001

**Our design fits everywhere**



- Automated photometric sample analysis for USP apparatus 1, 2, 5 and 6
- Unique Hollow Shaft™ sampling
- On-line UV-Vis and/or off-line sample collection into HPLC vials
- Automated media addition and pH change capabilities
- Autocompliance™ minimizes routine physical validation checks
- Controlled by 21 CFR Part 11 compliant WinSOTAX Advanced Dissolution Software

**The SOTAX AT 7*smart* on-line automated dissolution system with the Evolution 300 spectrophotometer from the Thermo Electron Corporation:**

- Long-life Xenon light source
- 8 and 16 position cell changers-control up to two dissolution systems with a single instrument
- Separate UV-Vis validation package available
- Automated performance and calibration tests

**Thermo**  
ELECTRON CORPORATION

**SOTAX**