

Test Certificate

Part: **Basket without Lock**
 Number: 02830
 Manufacturer: SOTAX
 Serial Number: 37001 - 37250

For SOTAX AT Dissolution Testing Systems



			USP/Ph.Eur. Tolerances	
Dimensions:	Basket rim O.D. Ø	25.0 mm	± 3.0	mm
	Height	37.0 mm	± 3.0	mm
	Wire Ø	0.25 mm	0.22 – 0.31	mm
	Wire openings	0.4 mm	0.36 - 0.44	mm
	Screen opening height	27.0 mm	± 1.0	mm
	Screen O.D. Ø	22.2 mm	± 1.0	mm
	Clear opening Ø	20.2 mm	± 1.0	mm
Material:	Stainless Steel			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
 Bruno Balli
 2013-11-29


 Signature