

Test Certificate

Part: **Basket without Lock**
 Number: 02830
 Serial Number: N130701 - N131050

For SOTAX AT Dissolution Testing Systems



| | | | USP/Ph.Eur. Tolerances | |
|-------------|-----------------------|---------|---------------------------|----|
| Dimensions: | Basket rim O.D. Ø | 25.0 mm | ± 3.0 | mm |
| | Height | 37.0 mm | ± 3.0 | mm |
| | Wire Ø | 0.25 mm | 0.22 – 0.31 | mm |
| | Wire openings | 0.4 mm | 0.36 - 0.44 | mm |
| | Screen opening height | 27.0 mm | ± 1.0 | mm |
| | Screen O.D. Ø | 22.2 mm | ± 1.0 | mm |
| | Clear opening Ø | 20.2 mm | ± 1.0 | mm |
| Material: | Stainless Steel | | | |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
 Bruno Balli
 2014-06-05

Signature