

# Test Certificate

Part: **Paddle with Suction Opening**  
 Number: 02985-2  
 Serial Number: 20751 - 21000

For SOTAX AT Dissolution Testing Systems



|             |        |              |                | USP/Ph.Eur.<br>Tolerances |
|-------------|--------|--------------|----------------|---------------------------|
| Dimensions: | Paddle | Height       | 19.0 mm        | ± 0.5 mm                  |
|             |        | Upper length | 74.5 mm        | ± 0.5 mm                  |
|             |        | Lower length | 42.0 mm        | ± 1.0 mm                  |
|             |        | Size         | 4.0 mm         | ± 1.0 mm                  |
|             |        | Radius       | 1.2 mm         | ± 0.2 mm                  |
|             |        | Radius       | 41.5 mm        | ± 1.0 mm                  |
|             | Shaft  | Diameter     | 9.4 to 10.1 mm | - -                       |

Material: **Stainless Steel (electrolytic polished)**

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
 Bruno Balli  
 2015-04-16

  
 Signature