


# Test Certificate

Part: **Vessel Mini Stirrer**  
 Number: 10583  
 Serial Number: 57763 - 57829

For SOTAX AT Dissolution Testing Systems



|             |   |           |          |                  |
|-------------|---|-----------|----------|------------------|
| Dimensions: | A | Outside Ø | 124.5 mm | + 0.3 / - 0.2 mm |
|             | B | Inside Ø  | 100.0 mm | ± 0.5 mm         |
|             | C | Outside Ø | 50.0 mm  | ± 1.0 mm         |
|             | D | Inside Ø  | 40.0 mm  | ± 1.0 mm         |
|             | E | Height    | 211.5 mm | ± 0.5 mm         |

With hemispherical bottom

Material: Borosilicate glass

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

SOTAX Quality Management  
 Bruno Balli  
 2014-10-24

  
 Signature