

# Test Certificate

Part: **Timer Flowability**

Number: 231532

For SOTAX PF1



Timer Accuracy: CUB7TV: 0.03% +100 m/sec per RUN terminal activation  
CUB7TC low freq/snk setup: 0.03% +1 msec per RUN terminal activation  
CUB7TC high freq/snk setup: 0.03% -1 msec per RUN terminal activation

Note: POWER: Non-replaceable internal 3.6 VDC lithium battery provides 7 years of typical continuous operation (high count speeds in SNK mode & extreme ambient temperatures will decrease battery life, use of SRC mode can extend battery life)

The above-mentioned equipment was checked of his accuracy. All relevant measurements show the requisite accuracy, and have passed the factory examination.

SOTAX Quality Management  
Bruno Balli  
2018-08-16

  
Signature