



Part

Paddle with Suction Opening

Part Number

02985-2

Serial Number

N13004838 – N13004933



Test Certificate

				USP / Ph.Eur. Tolerances
Dimensions	Paddle	Height	19.0 mm	± 0.5 mm
		Upper edge length	74.5 mm	± 0.5 mm
		Lower edge length	42.0 mm	± 1.0 mm
		Thickness	4.0 mm	± 1.0 mm
		Radius	1.2 mm	± 0.2 mm
		Radius	41.5 mm	± 1.0 mm
	Shaft	Diameter	9.4 – 10.1 mm	--
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2025-07-10

Buner Cornelius  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland