



Part Paddle with Stirring Shaft
Part Number 07525-04
Serial Number N16010931 - N16011074



Test Certificate

				USP / Ph.Eur. Tolerances
Dimensions	Paddle	Height	19.0 mm	± 0.5 mm
		Upper edge length	74.5 mm	± 0.5 mm
		Lower edge length	42.0 mm	± 1.0 mm
		Thickness	4.0 mm	± 1.0 mm
		Radius	1.2 mm	± 0.2 mm
		Radius	41.5 mm	± 1.0 mm
	Shaft	Diameter	9.4 – 10.1 mm	--
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2022-12-07

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland