



Part Basket Holder with Stirring Shaft  
Part Number 07558-06  
Serial Number N18005384 - N18006298



# Test Certificate

|                   |   |               | USP / Ph.Eur.<br>Tolerances |
|-------------------|---|---------------|-----------------------------|
| <b>Dimensions</b> | Shaft Ø                                 | 9.4 – 10.1 mm | - -                         |
|                   | Cap Ø                                   | 25.0 mm       | ± 3.0 mm                    |
|                   | Cap size                                | 5.1 mm        | ± 0.5 mm                    |
|                   | Clear opening                           | 20.2 mm       | ± 1.0 mm                    |
|                   | Vent hole Ø                             | 2.0 mm        | ± 0.5 mm                    |
| <b>Material</b>   | Stainless Steel (electrolytic polished) |               |                             |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2020-01-31

Bruno Balli  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland