



# Test Certificate

Part Suction Head with Stirring Shaft  
Part Number 07570  
Serial Number N15002308 - N15002409



				USP / Ph.Eur. Tolerances
<b>Dimensions</b>	Shaft	Diameter	9.4 – 10.1 mm	--
<b>Material</b>	Suction head PVC			
	Shaft Stainless Steel (electrolytic polished)			
	Suction tube PTFE			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2021-09-07

Bruno Balli  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland