




# Test Certificate

Part Vessel Mini Stirrer  
Part Number 10583  
Serial Number 84822 - 84835



		Tolerances	
Dimensions	Vessel	Outside Ø A	124.5 mm + 0.3 / - 0.2 mm
		Inside Ø B	100.0 mm ± 0.5 mm
		Outside Ø C	50.0 mm ± 1.0 mm
		Inside Ø D	40.0 mm ± 1.0 mm
		Height E	211.5 mm ± 0.5 mm
		With hemispherical bottom	
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

Aesch, 2021-01-20

Bruno Balli  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland