




Test Certificate

Part 1 Liter Vessel
Part Number 14143-01
Serial Number 103411 - 103921



		USP / Ph.Eur. Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm
		Height b	200 mm
		With hemispherical bottom	160 – 210 mm
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2025-08-14


06.11.2025 11:16

Michael Schmidlin
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland