




Test Certificate

PartAutomated Disintegration Basket
(for 6 Tubes)

Part Number23130-01

Serial NumberSK6.07860 – SK6.07869



| | | USP / Ph.Eur | Tolerances |
|------------|--------------------|--------------|------------|
| Dimensions | Bottom plate Ø | 90.0 mm | ± 2.0 mm |
| | Bottom thickness | 6.75 mm | ± 1.75 mm |
| | Distance to sieve | 6.75 mm | ± 1.75 mm |
| | 6 Holes Ø | 24 mm | ± 2.00 mm |
| Material | PMMA, POM, PVC, PA | | |

USP / Ph.Eur conformity (acc. USP701 / Ph Eur 2.9.1) :

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

The design of the basket-rack assembly may be varied somewhat, provided the specifications for the glass tubes and the screen mesh size are maintained.
(separate Certificates for Glass tube #12514 and Sieve plate # 23440-01)

Aesch, 2024-07-04

Michael Schmidlin
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland