




Part Automated Disintegration Disk
 Part Number 23444-02



**Ph.Eur.
Tolerances**

Dimensions	A	Disk outside Ø	31.40 mm	± 0.13 mm
	B	Thickness	15.30 mm	± 0.15 mm
	C	7 Holes Ø	3.15 mm	± 0.1 mm
	D	Internal bolt circle	8.40 mm	± 0.1 mm
Specific gravity	ϕ 1.18 – 1.20 kg/dm ³			--
Material	PMMA			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2021-03-10

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate