

Certificate of Compliance



Part SAM
Part Number 15832 (SR)
13135 (S)
13136 (M)
13137 (L)
Year of Construction from 2014-10

This document certifies that the above mentioned apparatus was assembled and successfully tested at the SOTAX AG facility in Aesch, Switzerland.

- DQ = Design qualification
- Visual check of system
- Electrical safety check
- Mechanical testing
- Calibration
- System testing
- Compatibility with the valid USP and Ph. Eur. requirements

We confirm that all parts in contact with the dissolution media are made of inert material such as:

Tubing	PTFE / PEEK
Tubing connector	PVDF
Valve body	FFKM
Syringe, Vials, Tubes	Glass
Needles & cannula	Stainless steel: 1.4401
Seals, Fittings, Gaskets	Viton (FPM)

Test data for this apparatus is available for inspection on request.

Aesch, 2019-08-05

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland