

Test Certificate


Part: **1 Liter Vessel with Connection for Drain Valve**

Number: 04871

Manufacturer: SOTAX

Serial Number: 54096 - 54137

For SOTAX AT 70*smart* Dissolution Testing Systems



USP/Ph.Eur.

Dimensions:	Inside Ø	a	100 mm	Tolerances	98 – 106 mm
	Height	b	200 mm		160 – 210 mm

With hemispherical bottom

Material: Borosilicate glass / Stainless Steel

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
Bruno Balli
2013-11-26


Signature