

Test Certificate



Part Paddle with Stirring Shaft
Part Number 07525-05
Serial Number P12 001221 – P12 001435



		USP / Ph.Eur. Tolerances		
Dimensions	Paddle	Height	19.0 mm	± 0.5 mm
		Upper edge length	74.5 mm	± 0.5 mm
		Lower edge length	42.0 mm	± 1.0 mm
		Thickness	4.0 mm	± 1.0 mm
		Radius	1.2 mm	± 0.2 mm
		Radius	41.5 mm	± 1.0 mm
	Shaft	Diameter	9.4 – 10.1 mm	--
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2023-08-31

31-08-2023 11:54:54

SOTAX Quality Management

C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India