

Test Certificate

Part: **Paddle with Stirring Shaft**
 Number: 07525
 Serial Number: 14001 – 14250

For SOTAX AT Dissolution Testing Systems



| | | | | |
|-------------|--------|--------------|----------------|---------------------------|
| | | | | USP/Ph.Eur. Tolerances |
| Dimensions: | Paddle | Height | 19.0 mm | ± 0.5 mm |
| | | Upper length | 74.5 mm | ± 0.5 mm |
| | | Lower length | 42.0 mm | ± 1.0 mm |
| | | Size | 4.0 mm | ± 1.0 mm |
| | | Radius | 41.5 mm | ± 1.0 mm |
| | Shaft | Diameter | 9.4 to 10.1 mm | - - |

Material: Stainless Steel (electrolytic polished)

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
 Bruno Balli
 2013-03-14


 Signature