



Part **Basket Holder with Stirring Shaft**  
Part Number **07558-07**  
Serial Number **P13 001931 – P13 002188**



			USP / Ph.Eur. Tolerances
<b>Dimensions</b>	Shaft Ø	9.4 – 10.1 mm	--
	Cap Ø	25.0 mm	± 3.0 mm
	Cap size	5.1 mm	± 0.5 mm
	Clear opening	20.2 mm	± 1.0 mm
	Vent hole Ø	2.0 mm	± 0.5 mm
<b>Material</b>	Stainless Steel (electrolytic polished)		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2024-08-08

Moreshwar Pulekar  
SOTAX Quality Management  
C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India

# Test Certificate