



Test Certificate

Part Suction Head with Stirring Shaft
 Part Number 07570
 Serial Number N15002410 - N15002675



			USP / Ph.Eur. Tolerances
Dimensions	Shaft	Diameter	9.4 – 10.1 mm - -
Material	Suction head PVC		
	Shaft Stainless Steel (electrolytic polished)		
	Suction tube PTFE		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2021-12-28

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland