

Test Certificate

Part: **Suction Head with Filter Insert (500 ml /900 ml)**
 Number: 07760
 Serial Number: 9701 - 9721

For SOTAX AT Dissolution Testing Systems



			USP/Ph.Eur. Tolerances
Dimensions:	Shaft Ø	9.4 to 10.1 mm	--
Material:	Suction head	PVC	
	Shaft	Stainless Steel	
	Suction tube	PTFE	

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
 Bruno Balli
 2013-11-01


 Signature