




Part 1 Liter Vessel
Part Number 14143-11
Serial Number P35 001305 - P35 001358



Test Certificate

		USP / Ph.Eur. Tolerances		
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	98 – 106 mm
		Height b	168.2 mm	160 – 210 mm
		With hemispherical bottom		
Material	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2022-02-14

14-02-2022 16:44:20

Moreshwar Pulekar
SOTAX Quality Management
C-64 | First Floor | TTC Industrial area | Turbhe | Navi Mumbai 400705 | India