




Part 1 Liter Vessel
 Part Number 14143-11
 Serial Number P35 001846 – P35 001958



Test Certificate

		USP / Ph.Eur. Tolerances		
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	98 – 106 mm
		Height b	168.2 mm	160 – 210 mm
		With hemispherical bottom		
Material	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Navi Mumbai, 2022-12-06

06-12-2022 15:49:26

Moreshwar Pulekar
 SOTAX Quality Management
 C-64 | First Floor | TTC Industrial area | Turbhe | Navi Mumbai 400705 | India