




Part 1 Liter Vessel  
 Part Number 14143-11  
 Serial Number P35 003318 – P35 003494



# Test Certificate

		USP / Ph.Eur. Tolerances	
<b>Dimensions</b>	Vessel	Inside Ø a	102.9 to 105.1 mm
		Height b	200 mm
		With hemispherical bottom	160 – 210 mm
<b>Material</b>	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2024-12-19

Moreshwar Pulekar  
 SOTAX Quality Management  
 C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India