




Part China Vessel 250 mL  
Part Number 14536-01  
Serial Number 88536 - 88722



|            |                    | PPRC Tolerances           |        |        |
|------------|--------------------|---------------------------|--------|--------|
| Dimensions | Vessel             | Inside Ø a                | 62 mm  | ± 3 mm |
|            |                    | Height b                  | 126 mm | ± 6 mm |
|            |                    | With hemispherical bottom |        |        |
| Material   | Borosilicate glass |                           |        |        |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with PPRC Pharmacopeia requirements (at date of approval, see below).

Aesch, 2022-03-11

Bruno Balli  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

# Test Certificate