




Part 1 Liter Vessel inactinic brown
Part Number 14537-02
Serial Number 98198 - 100116



			Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	--
		Height b	186.9 mm	± 0.7 mm
		With hemispherical bottom		
Material	Borosilicate glass			

Aesch, 2024-08-13

Tayyib Özüdogru
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate