

# Test Certificate

Part: **Basket without Lock**  
 Number: 02830  
 Manufacturer: SOTAX  
 Serial Number: 34251 - 34500

For SOTAX AT Dissolution Testing Systems



		USP/Ph.Eur. Tolerances	
Dimensions:	Basket rim O.D. Ø	25.0 mm	± 3.0 mm
	Height	37.0 mm	± 3.0 mm
	Wire Ø	0.28 mm	0.22 – 0.31 mm
	Wire openings	0.4 mm	0.36 - 0.44 mm
	Screen opening height	27.0 mm	± 1.0 mm
	Screen O.D. Ø	22.2 mm	± 1.0 mm
	Clear opening Ø	20.2 mm	± 1.0 mm
	Material:	Stainless Steel	

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
 Bruno Balli  
 2013-07-30

  
 Signature