

# Test Certificate

Part: **Suction Head with Stirring Shaft**  
 Number: 07570  
 Serial Number: N15000544 - N15000867

For SOTAX AT Dissolution Testing Systems



			USP/Ph.Eur. Tolerances
Dimensions:	Shaft Ø	9.4 to 10.1 mm	--
Material:	Suction head	PVC	
	Shaft	Stainless Steel	
	Suction tube	PTFE	

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
 Bruno Balli  
 2018-07-10

  
 Signature