

# Test Certificate

Part: **Suction Head with Filter Insert (500 ml /900 ml)**  
 Number: 07760  
 Serial Number: 9785 - 9839

For SOTAX AT Dissolution Testing Systems



|             |              |                 |                           |
|-------------|--------------|-----------------|---------------------------|
|             |              |                 | USP/Ph.Eur.<br>Tolerances |
| Dimensions: | Shaft Ø      | 9.4 to 10.1 mm  | --                        |
| Material:   | Suction head | PVC             |                           |
|             | Shaft        | Stainless Steel |                           |
|             | Suction tube | PTFE            |                           |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
 Bruno Balli  
 2015-03-10

  
 Signature