


# Test Certificate

Part: **Automated Disintegration Disk  
(for Capsules)**

Number: 08615

For SOTAX DT Tablet Disintegrators (6 Glass Basket)



Specific gravity

$\zeta 1.18 - 1.20 \text{ kg/dm}^3$

Material: PMMA transparent

All Dimensions in mm

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
Bruno Balli  
2017-09-25

  
Signature