


Test Certificate

Part: **Disintegration Disk**

Number: 08826

For SOTAX DT Tablet Disintegrators (3 Glass Basket)



Dimensions:	a	Disk outside Ø	31.4 mm	Ph.Eur. Tolerances ± 0.13 mm
	b	Thickness	15.3 mm	± 0.15 mm
	c	7 Holes Ø	3.15 mm	± 0.1 mm
	d	Internal bolt circle	8.4 mm	± 0.1 mm
		Specific gravity	ζ 1.18 – 1.20 kg/dm ³	--

Material: PMMA transparent

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with Ph. Eur. 2.9.1.-2 Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
Bruno Balli
2017-08-21


Signature