

Test Certificate

Part: **Protection Ring**
Number: 11148 and Set 11252 (7 pc)

For the Flange Protection of SOTAX Glass Vessels



Dimensions:	Outside Ø	125.6 mm	± 0.5 mm
	Inside Ø	124.0 mm	± 0.1 mm
	Inside Ø	40.0 mm	+ 0 / -0.3 mm
	Height	5.95 mm	± 0.1 mm

Max. temperature 60 °C

Material: PP

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

SOTAX Quality Management
Bruno Balli
2014-04-17


Signature