

# Certificate of Compliance


**SOTAX**

Solutions for Pharmaceutical Testing

SOTAX AG, Nordring 1, Aesch 4147, Switzerland

Part: **1 Liter Vessel**  
Number: 13411-01  
Serial Number: 72295 to 72315

For SOTAX AT Dissolution Testing Systems



Dimensions:	Inside Ø	a	100 mm	± 0.3 mm
	Height	b	186.9 mm	± 0.7 mm

With hemispherical bottom

Material: Borosilicate glass

SOTAX Quality Management  
Bruno Balli  
2017-10-05

  
Signature