


Test Certificate

Part: **Automated Disintegration Disk**

Number:	Variant	Dimension X
	18312-01	0.3 mm
	18312-02	0.5 mm
	18312-03	1.5 mm
	18312-04	2.0 mm
	18312-05	2.5 mm

Material: PMMA transparent
Specific gravity: ζ 1.18 – 1.20 kg/dm³

For SOTAX DT Tablet Disintegrators (6 Glass Basket)



All Dimensions in mm

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph.Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
Jürg Suter
2018-05-03


Signature