

Certificate of Compliance



Part FS Xtend™
Part Number 14460 (6-channels)
14470 (7-channels)
14480 (8-channels)
Year of Construction from 2013-08

This document certifies that the above mentioned apparatus was assembled and successfully tested at the SOTAX AG facility in Aesch, Switzerland.

- DQ = Design qualification
- Visual check of system
- Electrical safety check
- Mechanical testing
- Calibration
- System testing
- Compatibility with the valid USP and Ph. Eur. requirements

We confirm that all parts in contact with the dissolution media are made of inert material such as:

Tubes for media	PTFE
Filter holder top	PEEK
Filter holder below	PEEK
Seals	VITON

Test data for this apparatus is available for inspection on request.

Aesch, 2019-07-30

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland