



Part Suction Head with Stirring Shaft  
 Part Number 07570  
 Serial Number N15000001 - N15000336



# Test Certificate

			USP / Ph.Eur. Tolerances
<b>Dimensions</b>	Shaft	Diameter	9.4 – 10.1 mm
<b>Material</b>	Suction head PVC		
	Shaft Stainless Steel (electrolytic polished)		
	Suction tube PTFE		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2018-03-29

Bruno Balli  
 SOTAX Quality Management  
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland