# Test Certificate

| | |
|---|---|
| Part | Basket Holder with Stirring Shaft |
| Part Number | 07558-08 |
| Serial Number | N38000673 – N38000763 |



| | | | USP / Ph.Eur. Tolerances |
|---|---|---|---|
| **Dimensions** | Shaft Ø | 10 mm | 9.4 – 10.1 mm |
| | Cap Size | 4.85 mm | 5.1 ± 0.5 mm |
| | Clear opening | 20 mm | 20.2 ± 1.0 mm |
| | Vent hole Ø | 2 mm | 2 ± 1.0 mm |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2026-09-17

06.10.2026 07:22

Tayyib Özüdogru
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland



CT00000EN / 02