Test Certificate

Part	Basket Holder with Stirring Shaft
Part Number	07558-08
Serial Number	N38000764 – N38000784



			USP / Ph.Eur. Tolerances
Dimensions	Shaft Ø	10 mm	9.4 – 10.1 mm
	Cap Size	4.85 mm	5.1 ± 0.5 mm
	Clear opening	20 mm	20.2 ± 1.0 mm
	Vent hole Ø	2 mm	2 ± 1.0 mm

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2026-10-19

27.10.2026 08:46

Tayyib Özüdogru
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

CT00000EN / 02